**innoliving**  
easy tech, easy life

**ALCOOL TESTER**

**INN-991  
MANUALE D'USO**

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**ALCOHOL TESTER**

**INN-991  
USER MANUAL**

Thank you for purchasing the INNOLIVING alcohol tester, which uses a semiconductor gas sensor to detect the concentration of alcohol in the exhaled air and converts this measurement into a readable value that expresses the estimated concentration of alcohol in the blood.

**BEFORE USE, READ THIS USER MANUAL CAREFULLY AND KEEP IT FOR FUTURE REFERENCE.**

### WARNINGS

- The tester provides you with information to determine your blood alcohol level at the same time you use the test. Blood alcohol concentration changes as time passes and measurement may be affected by environmental conditions. This product is designed to be used as reference only in order to promote a responsible and

safe driving. It provides a purely indicative measurement, it doesn't replace or exclude possible tests conducted by competent authorities, which can have different readings. The tester readings have no legal value and do not relieve any individual responsibility. Improper storage of the product and/or improper use of the product may affect the accuracy of the measurement. You should not rely upon it as the sole basis to determine state of drunkenness or whether it is safe to drive a vehicle, operate equipment, or engage in dangerous activities. Innoliving Spa is not liable for the improper use of this product for any reason.

- Before use, make sure that the appliance is

intact without visible damage. If in doubt, do not use the product and contact the authorized after sales service

- After use, always store the product away from sources of heat, sunlight, humidity, sharp objects or similar.
- Do not use alcoholic substances such as paint or insecticides near the tester.
- Do not use chemicals to clean the tester. Never immerse the product in water or other liquids, the sensor could be damaged.
- Do not attempt to repair the product yourself, under penalty of invalidation of the warranty, but always contact the authorized after sales service.
- It takes about 20 minutes for alcohol

to be absorbed into the bloodstream by the digestive system; therefore, it is recommended to take the measurement 20 minutes after consuming alcoholic beverages. If the measurement is taken immediately after drinking, the results may reflect the concentration of alcohol in the drink and not the actual concentration in the blood

- To obtain more accurate measurement results, the interval between two consecutive measurements should be more than 3 minutes. Repeat the measurements 2 or 3 times to get a more accurate result.

**NOTE:** If the test is performed with the battery with a low charge level, the test value may differ from the real value. When the battery is

low, the device display will show "L" and the device will beep to signal the need to recharge the battery. When the appliance is charging, the indicator light will turn red; When fully charged the red light will turn off. On a full charge, the device can take about 200 measurements. The built-in battery is not replaceable. Do not attempt to remove it. The liquids that compose the battery are highly aggressive. In the event of any damage to the device resulting in the leakage of battery liquid, do not touch the product without adequate protection of the skin and eyes. Do not throw into fire, there is a risk of fire and explosion. Stop using the device immediately if you see swelling, leakage or other abnormalities.

### PRODUCT DESCRIPTION

### HOW TO USE THE ALCOHOL TESTER

The products work with a lithium built-in 120mAh 3,7V battery rechargeable through the USB cable supplied.

#### WARM-UP

1. Press the power button and hold for 2 seconds. You will hear a "beep" when turning on, the display is ON.
2. The word "WARM-UP" appears on the display, the device starts a 10-second countdown to pre-heating. Wait for pre-heating to complete.

#### NOTE:

Since the sensor surface will be oxidized when the test is not used for long time, when the tester is taken out of the package for the first time or the tester has not been used for a long time, the user needs to preheat the product as

many times as possible before performing the alcohol test.

### MEASUREMENT

1. At the end of the pre-heating, the device moves on to the measurement phase by blowing, the display shows the word "BLOWING" and a countdown from 10 to 0 begins.
2. Blow into the air intake during the countdown. Once the device detects the gas, a 3-second countdown begins and the audible alert will signal to stop blowing.
3. Read the result directly on the product display.
4. If the result of the detected alcohol concentration is 0.00%BAC(or 0.0g/l), it means that no alcohol has been detected and the LED indicator turns green.
5. If the result is higher than 0.01%BAC (0.1g/l) but less than 0.05%BAC (0.5g/l), the display will show "WARNING" and the LED indicator will turn yellow, accompanied by a beep.
6. If the result is higher than 0.05%BAC (0.5g/l), the display will show "DANGER" and the LED indicator will turn red, accompanied by a beep.

### SYMBOLS AND LED INDICATORS

**WARM UP:** after switching on the appliance, the word "WARM UP" appears on the display to indicate the pre-heating phase.

**BLOW:** when pre-heating is complete, the display indicates "BLOW" to start blowing on

the air intake.  
**RESULT:** when the measured value is displayed, the word "RESULT" and the corresponding measured value in g/l and %BAC at the same time appear in the display.  
**OFF:** "OFF" appears 10 seconds after the value is detected to signal auto-off.

Depending on the value detected, the light indicator takes on a different colour:

- Green: result lower than 0,1g/L
- Yellow: result from 0,1g/l to 0,4g/L
- Red: result higher than 0,5g/L

**CE** The product complies with all the applicable European regulations.

The battery used in this device must be disposed of in the special bins at the end of its life. The battery must be disassembled by the treatment center as it is not removable. The correct disposal of batteries helps preventing potentially negative consequences on the environment and human health.

**INFORMATION TO USERS** according to Legislative Decree No. 49 of March 14, 2014 " 2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE) " The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for

devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014

120mAh batteria ai polimeri di litio/Built-in Polymer Lithium ion Battery/  
Modello/Model:501520  
Fabbricante/Manufacturer: Shenzhen Juxinyuan Electronic Technology Co., Ltd Room 1804 Building, Zhongkegu Industrial Park, Zhonghuan Avenue, Shanzia Community, Pinghu Street, Longgang District Shenzhen China  
Importato da/Imported by:Innoliving Spa  
Via Merloni, 2/B 60131 Ancona Italy  
Data di produzione/Mfg date:  
Peso/weight:

### RATING LABEL

Innoliving Spa Via Merloni 2/B - 60131 Ancona - Italy  
**INN-991 ALCOHOL TESTER**  
5V== 1A BATTERIA INTEGRATA - BUILT-IN BATTERY  
3,7V== LITHIUM BATT. 120mAh  
MADE IN CHINA  
LOT n.

**innoliving**  
easy tech, easy life

Innoliving Spa  
Via Merloni, 2/B - 60131 Ancona Italy  
Tel 071.2133550 - [www.innoliving.it](http://www.innoliving.it)

MADE IN CHINA