INN
911

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IMPORTANT SAFETY RULES

- Read the instructions carefully before use and keep this manual for future reference.
- The appliance is intended for domestic use only.
- This product can be used by children over 8 years of age and by people with reduced physical, mental or sensory capabilities, or lack of experience and knowledge, only if under supervision or if they have received instructions regarding the use of the appliance in safely and understand the associated dangers.
- Make sure that children do not play with the appliance. Cleaning and maintenance operations must not be carried out by children without supervision.
- Never immerse the appliance or the base in water or other liquids.
- When the battery charger is connected to the mains, do not handle the appliance with wet hands or feet.
- Only use accessories recommended by the manufacturer.

- Do not use the appliance if the body or head is damaged.
- For hygienic reasons, if the product is used by more than one person, it is advisable to use a different head for each user.
- Do not pull the USB cable to disconnect the product from the power supply.
- Do not wrap the power cord around the appliance. Do not bend or twist the power cord.
- Disconnect the USB cable in the event of a power failure and before cleaning the product.
- Do not attempt to modify or repair the appliance yourself. Repairs must be carried out by an authorized service center or qualified personnel.
- This product is not a toy.
- Keep the product away from heat sources.
- The power cord cannot be replaced. If the cable is damaged the appliance must be scrapped.
- It is not an equipment for professional use in dental practices, use only to clean teeth

- and gums, in case of pain or bleeding contact your doctor.
- The battery included in the toothbrush cannot be removed by the user and the recovery and treatment center will remove it before disposal.
 - The toothbrush must be charged using only the charging base provided.

TECHNICAL SPECIFICATIONS

Battery Type: Lithium battery 3.7V/600mAh/2.22Wh

IP grade: IPX7

Spare toothbrush models: INN-90801new, INN-90802, INN-90803, INN-90804, INN-90805, INN-90806.

Charging base

LED indicator for charge level

CHARGER (not included):

Input: 100-240Vac, 50/60Hz, 1.4W

Output: DC 5V, 50mA

FUNCTIONS

- 3 cleaning modes
 - 1) Gentle (white light)
 - 2) Strong (green light)
 - 3) Pulse (blue light)
- Shutdown in 2 minutes
- Effective cleaning of lateral and medial teeth, gum spaces and crevices. Suitable for daily use for removing food residues in the oral cavity.

USAGE

- When the toothbrush is off, press the button to turn it on and hold to choose the operating mode. When idle, it will automatically turn off after 2 minutes. To manually turn it off while it is running, press the button.
- Insert the head into the toothbrush, moisten it and apply toothpaste. To avoid splashes, place

- the brush head behind the teeth before starting the movement of the bristles.
- Guide the toothbrush between your teeth to ensure effective cleaning of each one.
- Clean your teeth starting from the external surfaces, then the internal ones, ending with the chewing surfaces. Avoid applying excessive pressure when brushing.

CHARGING MODE

- Place the toothbrush on the charging base, making sure it is correctly connected. While charging, the red light will be on and will turn green when completed.
- During the charging process, the indicator light flashes and, once finished, remains lit.
- If the toothbrush is not used for a long period, the charge indicator may not appear immediately; wait 5-10 minutes for it to return to normal.
- When the battery is low, the red indicator flashes. Only use your toothbrush after charging.

CLEANING AND MAINTENANCE

After use, rinse the head under running water, remove excess water and store it in an upright position.

Do not immerse the toothbrush or its base in water. Remove the brush head and clean the main unit with a damp cloth.

It is recommended to replace the brush head every 3 months and every time you notice an excessive consumption of the bristles. Use only the heads suitable for this model.

After use store the product in a dry place away from heat sources.

HOW TO REPLACE THE BRUSH HEAD

1. Remove the head.
2. Insert the head directly into the holder as shown in the following figure.

RATING LABEL:

Read carefully the instruction manual.

IPX7

Degree of protection from contact with water.

Class III insulated electrical appliance

This product complies with all the applicable EU Directives.

INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 " 2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE) "The crossed-out dustbin symbol indicates that the product at the end of its life must be collected separately from other waste.

The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014

The user cannot remove the built-in battery, waste collection and treatment center shall extract and dispose of it. The batteries used in this device must be disposed of in the special bins at the end of their life. Please inform yourself about the local rules on separate collection of batteries.

The correct disposal of batteries helps preventing potentially negative consequences on the environment and human health.

innoliving
easy tech, easy life

Innoliving Spa
Via Merloni 2/B 60131 Ancona
Tel: 071 2133550
www.innoliving.it

MADE IN CHINA

Rev.04_09.2024